

implantswiss.com

NOVODENT SA
Y-Parc Technopole
Avenue des Sciences 11
CH-1400, Yverdon-les-Bains
Switzerland

01.2021 Rev.03

 **implant**swiss
SMART
CHOICE

IMPLANT TECHNOLOGY



» Content

Implantswiss Dental System	2
Swiss Master Precision	4
Product Traceability	5
Glass Package	6
Platform Color Coding	7
» Implantswiss System Specifications	8
S.R.A. Surface	12
Hybrid Design	13
Conical Connection with Hexagonal Index	14
Platform Switching	14
Self Cutting Edge	15
Hygienic Neck Design	15
» Implantswiss Implants	16
Bone Level Implant	18
Tissue Level Implant	19
» Surgery	20
Implantswiss Surgical Drills and Protocols	22
» Short Implant	28
Short Implant	29
Short Implant Surgical Kit	31
Short Implant Surgical Drills and Protocols	32
» Zygoma Implant	36
Zygoma Implant Surgical Kit	39
» Implantswiss Prosthetic Components	40
» Bone Level System Prosthetic Components	42
Bone Level Healing Abutments	45
Couple System	45
Solid System	47
Octa System	48
O-ring System	49
Locator System	50
Multi-Unit System	51
Bar System	52
» Tissue Level System Prosthetic Components	53
Tissue Level Healing Abutments	56
Couple System	56
Solid System	57
Octa System	58
O-ring System	59
Locator System	60
» Zygoma Dental Implant System Prosthetic Components	61
» Digital Solutions	62
Bone Level Digital Solutions	64
Tissue Level Digital Solutions	65
Multi Unit Abutment Digital Solutions	66
Octa Abutment Digital Solutions	66
Ti-Base Abutment	67
» Index Code	68

implantswiss

Switzerland

is celebrated and acknowledged as a world authority in

science, mechanical design and craftsmanship



As a country Switzerland has been one of the most important centers of inventiveness that guarantees

high precision and proficiency

Swiss products with impeccable mechanical designs and exquisite ingenuity can be used for a lifetime and add artistry to life.

Implantswiss Dental Implant System derives its **exceptional** structure from the inspirational and the magnificent Swiss artisan history.

Implantswiss is a strong, reliable and unique brand in dental Implantology due to its innovative design, surface technology and high quality.

Implantswiss is dedicated to Research and Development through collaborations with universities and it supports crucial research investment in any clinical environment. Implantswiss follows the latest trends and helps initiate extensive and contemporary projects which will lead into the future of implant dentistry.

Implantswiss Dental System


is the most prosperous and most versatile Implant system in Dentistry offering a range of solutions to all the simple and complex clinical situations.



SWISS MASTER PRECISION



Implantswiss products are precisely manufactured with state-of-the-art technology and diligent quality control mechanisms.




The perfectionist mindset of our Implantswiss team in each stage of the manufacturing and quality control process is the essence of our mission to deliver high quality and reliable products to the Dental world.

PRODUCT TRACEABILITY



Each Implantswiss product can be easily traced from production to delivery.

Each individual package has a barcode label with a lot number created by the ERP software system. All steps from production to stock, delivery and sales can be monitored with this barcode system; therefore, the traceability process is conducted flawlessly.



GLASS PACKAGE

GLASS PACKAGING

The preservation of nature and human health have always been the primary concern for Implantswiss. As a responsible supplier Implantswiss is aiming to limit the use of plastics in our packaging. Therefore, all Implantswiss Implants are packed in 100% natural and hygienic glass vials.

Since the glass structure is inert, glass packages are especially preferred over plastic containers for preserving the high value products.

High barrier property, ease of sterilization and nonreactivity to chemical substances, are few of the many advantages of glass packaging.

SMART PACKAGING

The implant can be simply removed from the glass bottle and placed in the osteotomy site using the specific fixture driver.

Each product package has a label containing all necessary information about the implant such as its type, length, width and platform color coding.

The platform color coding depicts the associated diameters of the implants and prosthetic components, the products are color coded as below.



PLATFORM COLOR CODING




TISSUE LEVEL

BONE LEVEL



» **IMPLANTSWISS
SYSTEM
SPECIFICATIONS**




» SRA* Surface


To ensure maximum Osseointegration after sandblasting, the surface of the implant is etched using Double Acid Etching technique.

Implantswiss S.R.A. surface has a unique homogeneous surface roughness that perfectly mimics the bone morphology.

* Sandblasted Rough Acid-etched Surface




3D-Area Parameters



Detail of the thread Bank and representation of the surface shape after removal of a 4th order polynomial

ISO 25178				
Height Parameters				
Sa	1,80	µm		
Sq	2,44	µm		
Feature Parameters				
S10z	15,1	µm	pruning = 5%	Ten point height
S5p	8,59	µm	pruning = 5%	Five point peak height
S5v	6,49	µm	pruning = 5%	Five point of height




Detail of the thread valley and representation of the surface shape after removal of a 4th order polynomial

ISO 25178				
Height Parameters				
Sa	1,25	µm		
Sq	1,60	µm		
Feature Parameters				
S10z	9,16	µm	pruning = 5%	Ten point height
S5p	3,39	µm	pruning = 5%	Five point peak height
S5v	5,76	µm	pruning = 5%	Five point of height

» Tapered form + Straight form =

Hybrid Design

It is a unique & innovative hybrid implant design with cylindrical (straight) configuration at the neck area and conical (tapered) configuration at the apex of the implant.



The top 3 mm of Implantswiss implant is straight downward from the neck which constitutes the straight form, however the remaining apical portion of implant constitutes the tapered form at a certain angle. The optimum combination of these two specifications creates our hybrid design.


The straight portion at the top of our unique hybrid design exerts minimal stress on the crucial crestal cortical bone leading to an optimal load distribution. The tapered portion of the design facilitates easy insertion of the implant at osteotomy site

» Conical connection with hexagonal index

creates a secure and steady implant-abutment connection

The universal conical connection and hexagonal index of our implant creates a firm and sound implant-abutment interface, which prevents the formation of micro gap and hence micro leakage.

In Implantswiss Dental Implant System, all abutments are compatible with all fixtures of any diameter (except for 3.3mm diameter) making it one of the most user-friendly implant systems for dentists.



» Platform Switching

Maintains the crestal alveolar bone around the neck area of the implant.

Micro-gap can be formed at the implant-abutment interface in Bone Level implants especially if they have a butt joint connection. Pathogens and their by-products in these areas can cause serious problems such as bone resorption and peri-implantitis.

Although Implantswiss Implants have a sealed conical connection, Platform switching further protects & maintains cervical bone by moving the implant-abutment interface away from the alveolar bone. Implantswiss Implants with Platform switch characteristic are safer against bacteria and microleakage, thus have a longer life span.



» Self Cutting Edge

Unique thread design that does not transfer stress to the bone.


Implantswiss implants have a particular thread design which does not transfer stress to the bone during insertion. Semi-aggressive thread design with cutting edge and self-cutting features provide ideal primary stability in any bone type without transferring unfavorable stress to the bone.



» Hygienic Neck Design

Easy to maintain the hygiene

Harmful pathogens can easily colonize in the micro threads at the cervical area of some implant designs. The absence of micro threads at the cervical area of Implantswiss Implants make it easier to maintain hygiene and plaque control. Due to this feature, plaque around the neck of the Implants is less retentive and hence less chances of bone resorption.



» IMPLANTSWISS
IMPLANTS




BONE LEVEL IMPLANT



Bone Level implants are designed to provide excellent outcomes in any bone quality and at any site.

Due to their unique design, Implantswiss bone level implants allow you to design an ideal smile in the aesthetic zone.




System » Bone Level
Material » Ti Grade 4
Driver » Bone Level driver,
 Mini Bone Driver for 33 implant
Torque Value » 35 Ncm

TISSUE LEVEL IMPLANT



Excellent soft and hard tissue preservation can be achieved successfully with the smooth neck, at the implant-abutment interface of our Tissue Level implant.


The implant is designed to provide excellent outcomes in any situation where the aesthetic outcome is not the primary concern.



System » Tissue Level
Material » Ti Grade 4
Driver » Tissue level driver
Torque Value » 35 Ncm

» SURGERY

BONE LEVEL - TISSUE LEVEL





IMPLANTSWISS FIXTURE DRILLS & PROTOCOLS

SURGICAL DRILLS & PROTOCOLS

Self-stoppered surgical drills eliminate the risk of complications during the preparation of the osteotomy site and insertion of the implant. Implantswiss drills are recommended to be used at 200-300 rpm to prevent overheating and reduce mechanical trauma, therefore shortening the osseointegration time of the implant. The drills are designed to harvest the autogenous bone when used below 50RPM without irrigation, which can be used as a graft material in regenerative procedures.

The Implantswiss Surgical Kit has color-coded drills for each implant type available in the system.


The Implantswiss Surgical Kit has various drills with different diameters and lengths. These color-coded drills are designed to simplify the surgical procedures, resulting in shorter surgery time and less stress for your patients.




**ø3.7
Fixture**




D1 Bone



D2-3 Bone




D4 Bone




**ø4.3
Fixture**




D1 Bone



D2-3 Bone




D4 Bone




**ø4.8
Fixture**




D1 Bone



D2-3 Bone




D4 Bone




**ø5.5
Fixture**




D1 Bone



D2-3 Bone



D4 Bone



» **IMPLANTSWISS
SHORT IMPLANTS**



SHORT IMPLANT

Short implants are considered as an alternative treatment to advanced surgical interventions in case of insufficient bone height, either above the mandibular canal or below the maxillary sinuses and aim to regenerate hard tissues.

Due to the aggressive thread design, Bone level short implants can preferably be used in the maxilla where the alveolar bone is relatively soft, whereas Tissue level short implants can be a good choice for the mandible.

BONE LEVEL

Length	Ø4.3	Ø4.8	Ø5.5
5.0 mm	S-BMFSR4305	S-BMFSR4805	S-BMFSR5505
6.0 mm	S-BMFSR4306	S-BMFSR4806	S-BMFSR5506


TISSUE LEVEL

Length	Ø4.3	Ø4.8	Ø5.5
4.0 mm	S-TFSR4304	S-TFSR4804	S-TWFSR5504
5.0 mm	S-TFSR4305	S-TFSR4805	S-TWFSR5505
6.0 mm	S-TFSR4306	S-TFSR4806	S-TWFSR5506

GH > 265 mm


SHORT IMPLANT SURGICAL KIT

The short surgical kit is recommended for both Bone and Tissue Level short implants.



SHORT IMPLANT DRILLS & PROTOCOLS

The Implantswiss Short Implant Surgical Kit has a specific drill for each implant available in the system.




Implantswiss Short Implant Surgical Kit has various color-coded drills with different diameters and lengths. These drills are designed to simplify the surgical procedures resulting in shorter surgery time and less stress for the patients.




BONE LEVEL

D4 / D3 Bone



TISSUE LEVEL

D2 / D1 Bone




Bone Level Short Implants are recommended to be used only in D3/D4 and Tissue Level Implants only in D1/D2 bone.

**ø4.8
Fixture**



BONE LEVEL

D4 / D3 Bone




**ø5.5
Fixture**




BONE LEVEL

D4 / D3 Bone




TISSUE LEVEL

D2 / D1 Bone



TISSUE LEVEL

D2 / D1 Bone



» **IMPLANTSWISS
ZYGOMA IMPLANTS**



ZYGOMA IMPLANT

Zygomatic implants are designed for complex cases where conventional implants cannot be used and for patients who want to benefit from immediate loading protocols. Even in the presence of minimum residual maxillary bone height, zygomatic implants provide excellent primary stability.

Implantswiss Zygomatic Implant offers three different product options.

ZYGOMA IMPLANT

ø4.3



35 mm	S-ZYG35
37.5 mm	S-ZYG37.5
40 mm	S-ZYG40
42.5 mm	S-ZYG42.5
45 mm	S-ZYG45
47.5 mm	S-ZYG47.5
50 mm	S-ZYG50
52.5 mm	S-ZYG52.5
55 mm	S-ZYG55
57.5 mm	S-ZYG57.5
60 mm	S-ZYG60

ATROPHIC MAXILLARY ZYGOMA IMPLANT

ø4.3



35 mm	S-ZAM35
37.5 mm	S-ZAM37.5
40 mm	S-ZAM40
42.5 mm	S-ZAM42.5
45 mm	S-ZAM45
47.5 mm	S-ZAM47.5
50 mm	S-ZAM50
52.5 mm	S-ZAM52.5
55 mm	S-ZAM55
57.5 mm	S-ZAM57.5
60 mm	S-ZAM60


ONCOLOGICAL ZYGOMA IMPLANT

ø4.3



25 mm	SZ-ONC25
27.5 mm	SZ-ONC27.5
30 mm	SZ-ONC30
32.5 mm	SZ-ONC32.5

ZYGOMA IMPLANT SURGICAL KIT



SURGICAL DRILLS

Point drill Pilot drill 3.5 Pilot drill Long Twist drill 2.5 Short Twist drill 2.5 Long Twist drill 2.9 Short Twist drill 2.9 Long Twist drill 3.5 Short Twist drill 3.5 Long

















» IMPLANTSWISS
PROSTHETIC COMPONENTS




















BONE LEVEL















COUPLE SYSTEM

FLOW	FIXTURE	HEALING	IMPRESSION COPING	LAB ANALOG	ABUTMENT	SCREW
 		 Closing Screw	 Pick-up Hex	 LAB ANALOG	 Couple Hex	 SCREW
		 Healing Abutment	 Transfer Hex		 Angled 15° 25°	
TOOL	 Ratchet Driver	 Hex Screw Driver	 Hex Screw Driver			 Hex Screw Driver













OCTA SYSTEM

FLOW	FIXTURE	HEALING	SUB OCTA ABUTMENT	IMPRESSION COPING	LAB ANALOG	LAB CYLINDER	SCREW
 		 Closing Screw	 Octa Healing Cap	 Pick-up	 LAB ANALOG	 Octa	 SCREW
		 Healing Abutment	 Octa Abutment	 Transfer		 N-Octa	
TOOL	 Ratchet Driver	 Hex Screw Driver	 Octa Driver			 Octa/N-Octa	 Hex Screw Driver

SOLID SYSTEM

FLOW	FIXTURE	HEALING	ABUTMENT	IMPRESSION CAP	LAB ANALOG	PLASTIC COPING
 		 Closing Screw	 Solid	 IMPRESSION CAP	 LAB ANALOG	 Single
		 Healing Abutment	 Cap			 Bridge
TOOL	 Ratchet Driver	 Hex Screw Driver	 Solid Driver			 Hex Screw Driver

O-RING SYSTEM

FLOW	FIXTURE	HEALING	O-RING ABUTMENT	O-RING LAB ANALOG	RETAINER
 		 Closing Screw	 O-RING ABUTMENT	 O-RING LAB ANALOG	 S-ORR
		 Healing Abutment			 S-ORRc
TOOL	 Ratchet Driver	 Hex Screw Driver	 O-Ring Driver		 S-ORRo

LOCATOR SYSTEM

FLOW	FIXTURE	HEALING	LOCATOR ABUTMENT	LOCATOR LAB ANALOG	INSERTER
TOOL	 Ratchet Driver	 Hex Screw Driver	 Locator Driver		 S-LRI

MULTI-UNIT SYSTEM

FLOW	FIXTURE	HEALING	ABUTMENT	IMPRESSION COPING	LAB ANALOG	LAB CYLINDER
TOOL	 Ratchet Driver	 Hex Screw Driver	 Hex Screw Driver Multi-Unit Driver	 Hex Screw Driver		 Hex Screw Driver

BONE LEVEL HEALING ABUTMENTS

BONE LEVEL HEALING ABUTMENTS

	 ø3.7 <small>only for ø3.3 implant</small>	 ø4.5	 ø5.5	 ø6.5
1 mm	-	S-BHA4510	S-BHA5510	S-BHA6510
2 mm	-	S-BHA4520	S-BHA5520	S-BHA6520
2.5 mm	S-BMHA3725	-	-	-
3 mm	-	S-BHA4530	S-BHA5530	S-BHA6530
3.5 mm	S-BMHA3735	-	-	-
4 mm	-	S-BHA4540	S-BHA5540	S-BHA6540
5 mm	S-BMHA3750	S-BHA4550	S-BHA5550	S-BHA6550
7 mm	S-BMHA3770	-	-	-

System » Bone Level
Material » Ti Grade 5 ELI
Driver » Hex Screw Driver

BONE LEVEL SHORT IMPLANT HEALING ABUTMENTS

	 ø5.0	 ø6.0	 ø8.0
2 mm	S-BMHA5020	S-BMHA6020	S-BMHA8020
3 mm	S-BMHA5030	S-BMHA6030	S-BMHA8030

System » Bone Level
Material » Ti Grade 5 ELI
Driver » Hex Screw Driver




BONE LEVEL COUPLE SYSTEM

BONE LEVEL COUPLE ABUTMENTS

	 ø3.7 <small>only for ø3.3 implant</small>	 ø4.5	 ø5.5	 ø6.5
1 mm	S-BMCA37107	S-BCA45107	S-BCA55107	S-BCA65107
2 mm	S-BMCA37207	S-BCA45207	S-BCA55207	S-BCA65207
3 mm	S-BMCA37307	S-BCA45307	S-BCA55307	S-BCA65307
4 mm	S-BMCA37407	S-BCA45407	S-BCA55407	S-BCA65407
5 mm	S-BMCA37507	S-BCA45507	S-BCA55507	S-BCA65507





System » Bone Level
Material » Ti Grade 5 ELI
Torque Value » 25 Ncm
Driver » Hex Screw Driver

BONE LEVEL SHORT IMPLANT COUPLE ABUTMENTS

			
	Ø5.0	Ø6.0	Ø8.0
2 mm	S-BMCA5020	S-BMCA6020	S-BMCA8020
3 mm	S-BMCA5030	S-BMCA6030	S-BMCA8030




System » Bone Level
Material » Ti Grade 5 ELI
Torque Value » 25 Ncm
Driver » Hex Screw Driver

BONE LEVEL ANGLED ABUTMENT

				
	Ø3.7 <small>only for ø3.3 implant</small>	Ø4.5 <small>15° 25°</small>	Ø5.5 <small>15° 25°</small>	Ø6.5 <small>15° 25°</small>
1 mm (15°)	S-BMAA37151	2 mm 15° S-BAA45152	S-BAA55152	S-BAA65152
2 mm (15°)	S-BMAA37152	2 mm 25° S-BAA45252	S-BAA55252	S-BAA65252
4 mm (15°)	S-BMAA37154	4 mm 15° S-BAA45154	S-BAA55154	S-BAA65154
		4 mm 25° S-BAA45254	S-BAA55254	S-BAA65254





System » Bone Level // **Material** » Ti Grade 5 ELI // **Torque Value** » 25 Ncm // **Driver** » Hex Screw Driver

TEMPORARY ABUTMENTS

		
Ø3.7 <small>only for ø3.3 implant</small>	Ø4.5	Ø5.5
S-BMTA37	S-BTA45	S-BTA55



System » Bone Level
Material » PEEK
Torque Value » 20 Ncm
Driver » Hex Screw Driver

BONE LEVEL IMPRESSION COPING PICK UP

			
Ø3.7 <small>only for ø3.3 implant</small>	Ø4.5	Ø5.5	Ø6.5
S-BMICP37L S-BMICP37S	S-BICP45L S-BICP45S	S-BICP55L S-BICP55S	S-BICP65L S-BICP65S



System » Bone Level
Material » Ti Grade 5 ELI

BONE LEVEL SHORT IMPLANT IMPRESSION COPING PICK UP

	
Ø5.0	Ø6.0
S-BMSICP50	S-BMSICP60

System » Bone Level
Material » Ti Grade 5 ELI

BONE LEVEL IMPRESSION COPING TRANSFER

	
Ø3.7 <small>only for ø3.3 implant</small>	Ø4.5
S-BMICT37L S-BMICT37S	S-BIP45



System » Bone Level
Material » Ti Grade 5 ELI

BONE LEVEL SHORT IMPLANT IMPRESSION COPING TRANSFER

	
S-BMSICT50	S-BMSICT60

System » Bone Level
Material » Ti Grade 5 ELI




BONE LEVEL FIXTURE LAB ANALOG

	
S-BMFLA <small>only for ø3.3 implant</small>	S-BFLA

System » Bone Level
Material » Ti Grade 5 ELI


BONE LEVEL SOLID SYSTEM

BONE LEVEL SOLID ABUTMENTS

			
	Ø3.7 <small>only for ø3.3 implant</small>	Ø4.5	Ø5.5
1 mm	S-BMSA3710	S-BSA4510	S-BSA5510
2 mm	S-BMSA3720	S-BSA4520	S-BSA5520
3 mm	S-BMSA3730	S-BSA4530	S-BSA5530
4 mm	S-BMSA3740	S-BSA4540	S-BSA5540
5 mm	S-BMSA3750	S-BSA4550	S-BSA5550


System » Bone Level
Material » Ti Grade 5 ELI
Torque Value » 25 Ncm
Driver » Hex Screw Driver

BONE LEVEL SOLID ABUTMENT CAP




System » Bone Level Solid System
Material » POM

BONE LEVEL SOLID ABUTMENT IMPRESSION CAP



System » Bone Level Solid System
Material » Plastic

BONE LEVEL SOLID LAB ANALOG



System » Bone Level Solid System
Material » Ti Grade 5 ELI


BONE LEVEL OCTA SYSTEM

BONE LEVEL OCTA ABUTMENTS



System » Bone Level
Material » Ti Grade 5 ELI
Torque Value » : Abutment Screw
 25 Ncm Oclusal Screw 20 Ncm
Driver » Octa Driver

BONE LEVEL OCTA ABUTMENT HEALING CAP



System » Octa System
Material » Ti Grade 4

BONE LEVEL OCTA SYSTEM PLASTIC CYLINDER



System » Octa System
Material » POM

BONE LEVEL OCTA SYSTEM TITANIUM CYLINDER



System » Octa System
Material » Ti Grade 5 ELI

BONE LEVEL OCTA SYSTEM IMPRESSION COPING PICK UP



BONE LEVEL OCTA SYSTEM IMPRESSION COPING TRANSFER



System » Octa System
Material » Ti Grade 5 ELI


BONE LEVEL OCTA ABUTMENT LAB ANALOG



System » Octa System
Material » Ti Grade 5 ELI

BONE LEVEL O-RING SYSTEM

BONE LEVEL O-RING ABUTMENTS



System » Bone Level
Material » Ti Grade 5 ELI
Torque Value » 25 Ncm
Driver » O-Ring Driver

DIRECTIONAL RINGS



S-ORRa

O-RING LAB ANALOG



S-ORLA

System » O-Ring System
Material » Ti Grade 5 ELI

RETAINER



S-ORR S-ORRc S-ORG-Pcap S-ORG-Ycap S-ORG-Gcap

PINK Soft 900 Gr
YELLOW ExtraSoft 500 Gr
GREEN Elastic 350 Gr

BONE LEVEL LOCATOR SYSTEM

BONE LEVEL MINI LOCATOR ABUTMENTS

only for ø3.3 implant

1.5 mm	S-BMLA15
2.5 mm	S-BMLA25
3.5 mm	S-BMLA35
4.5 mm	S-BMLA45
5.5 mm	S-BMLA55
6.5 mm	S-BMLA65

BONE LEVEL LOCATOR ABUTMENTS

1.5 mm	S-BLA15
2.5 mm	S-BLA25
3.5 mm	S-BLA35
4.5 mm	S-BLA45
5.5 mm	S-BLA55
6.5 mm	S-BLA65



System » Bone Level
Material » Ti Grade 5 ELI
Torque Value » 25 Ncm
Driver » Locator Driver

BONE LEVEL LOCATOR LAB ANALOG



S-BLLA

LOCATOR MATRIX



S-LMht




BLACK CAP
Processing




VIOLET CAP
Rigid Retention 2.7 Kg



PINK CAP
Soft Retention 1.2 Kg



CLEAR CAP
Standard Retention 1.8 Kg



YELLOW CAP
Extra Soft Retention 0.6 Kg

S-LRI

BONE LEVEL MULTI UNIT SYSTEM

MULTI UNIT ABUTMENTS



ø4.8

only for ø3.3 implant

0.5 mm	S-BMMUA4805
1.5 mm	S-BMMUA4815
2.5 mm	S-BMMUA4825
3.5 mm	S-BMMUA4835
4.5 mm	S-BMMUA4845



ø4.8

0.5 mm	S-BMUA4805
1.5 mm	S-BMUA4815
2.5 mm	S-BMUA4825
3.5 mm	S-BMUA4835
4.5 mm	S-BMUA4845

System » Bone Level
Material » Ti Grade 5 ELI
Torque Value » : Abutment Screw
25 Ncm Oclusal Screw 20 Ncm
Driver » Multi Unit Driver

MULTI UNIT ANGLED ABUTMENTS



ø4.8 (17°)

1.0 mm	S-BMUA48171
2.5 mm	S-BMUA48172
3.5 mm	S-BMUA48173



ø4.8 (35°)

1.0 mm	S-BMUA48351
2.5 mm	S-BMUA48352
3.5 mm	S-BMUA48353

System » Bone Level
Material » Ti Grade 5 ELI
Torque Value » : Abutment Screw
25 Ncm Oclusal Screw 20 Ncm
Driver » 1,2 Hex Driver

MULTI UNIT ABUTMENT HEALING CAP



ø4.8

4 mm	S-MUAHCap
6 mm	S-MUAHCapL

System » Multi Unit System
Material » Ti Grade 5 ELI

MULTI UNIT IMPRESSION COPING PICK UP



ø4.8

S-BMUAICP

System » Multi Unit System
Material » Ti Grade 5 ELI

MULTI UNIT ABUTMENT PLASTIC CYLINDER



S-MUAPICyN

System » Multi Unit System
Material » POM

MULTI UNIT ABUTMENT TITANIUM CYLINDER



S-MUATICyN

System » Multi Unit System
Material » Ti Grade 5 ELI

MULTI UNIT LAB ANALOG



S-MUALA
Ø 4.8×14mm

System » Bone Level
Material » Ti Grade 5 ELI

BAR SYSTEM

BAR BURNOUT

S-IB



BAR CLIP YELLOW

S-IBClip - y
Medium Retention 1.8 Kg



TOOL FOR INSERTING CLIP



S-BMIA

BAR POSITIONING

S-IBP



BAR CLIP PINK

S-IBClip - p
Soft Retention 1.2 Kg



KEY FOR PARALLELOMETER



S-TFIC

BAR HOUSING

S-IBH



BAR CLIP RED

S-IBClip - r
Extra Soft Retention 0.8 Kg



TISSUE LEVEL

COUPLE SYSTEM

FLOW	FIXTURE	HEALING	ABUTMENT	LAB ANALOG	ABUTMENT	SCREW
 	 Healing Abutment	 Cover Screw	 Pick-up Transfer		 Couple Angled Synocta	 Hex Screw Driver
TOOL	 Ratchet Driver	 Hex Screw Driver	 Hex Screw Driver			 Hex Screw Driver

SOLID SYSTEM

FLOW	FIXTURE	HEALING	ABUTMENT	IMPRESSION CAP	LAB ANALOG	PLASTIC COPING
 	 Healing Abutment	 Closing Screw	 Cap		 Lab Analog	 Crown Bridge
TOOL	 Ratchet Driver	 Hex Screw Driver	 Solid Driver Hex Screw Driver			 Hex Screw Driver

OCTA SYSTEM

FLOW	FIXTURE	HEALING	OCTA ABUTMENT	IMPRESSION COPING	LAB ANALOG	LAB CYLINDER	SCREW
 ø3.7 ø4.3-L06 ø4.3	 ø4.8-L06 ø4.8 ø5.5	 Cover Screw	 Octa Abutment	 Pick-up	 Lab Analog	 Octa N-Octa	 Screw
TOOL	 Ratchet Driver	 Hex Screw Driver	 Octa Driver			 Octa/N-Octa	 Hex Screw Driver

LOCATOR SYSTEM

FLOW	FIXTURE	HEALING	LOCATOR ABUTMENT	LOCATOR LAB ANALOG	INSERTER
 ø3.7 ø4.3-L06 ø4.3	 ø4.8-L06 ø4.8 ø5.5	 Closing Screw	 Locator Abutment	 Locator Lab Analog	 S-LMHt
TOOL	 Ratchet Driver	 Hex Screw Driver	 Locator Driver		 S-LRI

O-RING SYSTEM

FLOW	FIXTURE	HEALING	O-RING ABUTMENT	O-RING LAB ANALOG	RETAINER
 ø3.7 ø4.3-L06 ø4.3	 ø4.8-L06 ø4.8 ø5.5	 Closing Screw	 O-Ring Abutment	 O-Ring Lab Analog	 S-ORR S-ORRc
TOOL	 Ratchet Driver	 Hex Screw Driver	 O-Ring Driver		 S-ORR S-ORRc

TISSUE LEVEL HEALING ABUTMENTS

TISSUE LEVEL HEALING ABUTMENTS



3.5 mm S-THA35 S-THAW35
 5.5 mm S-THA55 S-THAW55

only for ø5.5 implant

System » Tissue Level
 Material » Ti Grade 4
 Torque Value » 20 Ncm
 Driver » Hex Screw Driver

TISSUE LEVEL COUPLE SYSTEM

TISSUE LEVEL COUPLE ABUTMENTS



1 mm S-TCA10 S-TCAW10
 2 mm S-TCA20 S-TCAW20
 3 mm S-TCA30 S-TCAW30
 4 mm S-TCA40 S-TCAW40

only for ø5.5 implant

System » Tissue Level
 Material » Ti Grade 4
 Torque Value » 25 Ncm
 Driver » Hex Screw Driver

TISSUE LEVEL ANGLED ABUTMENT



15° S-TAA15
 25° S-TAA25

System » Tissue Level
 Material » Ti Grade 4
 Torque Value » 25 Ncm
 Driver » Hex Screw Driver

TISSUE LEVEL SYNOCTA ABUTMENT



6 mm S-TSOA

System » Tissue Level
 Material » Ti Grade 4
 Torque Value » 25 Ncm
 Driver » Hex Screw Driver

TISSUE LEVEL IMPRESSION COPING PICK UP



OCTA S-TICP S-TICPW

only for ø5.5 implant

System » Tissue Level
 Material » Ti Grade 5 ELI

TISSUE LEVEL IMPRESSION COPING TRANSFER



OCTA S-TICT

System » Tissue Level
 Material » Ti Grade 5 ELI

TISSUE LEVEL FIXTURE LAB ANALOG



S-TFLA



S-TFLLA

only for ø5.5 implant

System » Tissue Level
 Material » Ti Grade 5 ELI

TISSUE LEVEL SOLID SYSTEM

TISSUE LEVEL SOLID ABUTMENTS



ø3.5

4 mm S-TSA40
 5 mm -
 5.5 mm S-TSA55
 7 mm S-TSA70



ø4.5

only for ø5.5 implant

S-TSAL40
 S-TSAL50
 -
 S-TSAL70

System » Tissue Level
 Material » Ti Grade 4
 Torque Value » 25 Ncm
 Driver » Solid Driver

TISSUE LEVEL SOLID ABUTMENT POSITIONING CYLINDER & IMPRESSION CAP



S-TSAPoCy



S-TSAICap

System » Tissue Level Solid System
 Material » Plastic



TISSUE LEVEL SOLID PLASTIC COPING



CROWN S-TSAPICc
 BRIDGE S-TSAPICb

System » Tissue Level Solid System
 Material » POM

TISSUE LEVEL SOLID ABUTMENT LAB ANALOG

		
	only for ø5.5 implant	
4 mm	S-TSALA40	S-TSALLA40
5.5 mm	S-TSALA55	S-TSALLA55
7 mm	S-TSALA70	S-TSALLA70

System » Tissue Level Solid System
Material » Ti Grade 5 ELI

TISSUE LEVEL OCTA SYSTEM

TISSUE LEVEL OCTA ABUTMENT



System » Tissue Level
Material » Ti Grade 4
Torque Value » : Abutment Screw 25 Ncm Oclusal Screw 20 Ncm
Driver » Octa Driver

TISSUE LEVEL OCTA ABUTMENT HEALING CAP

	
4 mm	S-OAHCap
6 mm	S-OAHCapL

System » Octa System
Material » Ti Grade 5 ELI

TISSUE LEVEL OCTA ABUTMENT PLASTIC CYLINDER OCTA / NON-OCTA

		
OCTA	S-OAPICy	
NON-OCTA	S-OAPICyN	

System » Octa System
Material » POM

TISSUE LEVEL OCTA ABUTMENT TITANIUM CYLINDER OCTA / NON-OCTA

	
OCTA	S-OATICy
NON-OCTA	S-OATICyN

System » Octa System
Material » Ti Grade 5 ELI

TISSUE LEVEL OCTA SYSTEM IMPRESSION COPING PICK UP OCTA / NON-OCTA

	
	ø5.5
OCTA	S-ICOAP
N-OCTA	S-ICOAPN

System » Octa System
Material » Ti Grade 5 ELI

TISSUE LEVEL OCTA SYSTEM IMPRESSION COPING TRANSFER OCTA / NON-OCTA

	
	ø5.0
OCTA	S-ICOAT
N-OCTA	S-ICOATN

System » Octa System
Material » Ti Grade 5 ELI


TISSUE LEVEL OCTA LAB ANALOG



System » Octa System
Material » Ti Grade 5 ELI

TISSUE LEVEL O-RING SYSTEM

TISSUE LEVEL O-RING ABUTMENTS

	
0.5 mm	S-TORA05
2 mm	S-TORA20
4 mm	S-TORA40

System » Tissue Level
Material » Ti Grade 4
Torque Value » 25 Ncm
Driver » O-Ring Driver

DIRECTIONAL RINGS



O-RING LAB ANALOG



System » O-Ring System
Material » Ti Grade 5 ELI

RETAINER



PINK
Soft 900 Gr

YELLOW
ExtraSoft 500 Gr

GREEN
Elastic 350 Gr

S-ORR

S-ORRc

S-ORG-Pcap

S-ORG-Ycap

S-ORG-Gcap

TISSUE LEVEL LOCATOR SYSTEM

TISSUE LEVEL LOCATOR ABUTMENTS



1.5 mm	S-TLA15
2.5 mm	S-TLA25
3.5 mm	S-TLA35
4.5 mm	S-TLA45
5.5 mm	S-TLA55
6.5 mm	S-TLA65

System » Tissue Level
Material » Ti Grade 5 ELI
Torque Value » 25 Ncm
Driver » Locator Driver

TISSUE LEVEL LOCATOR LAB ANALOG



S-TLLA

LOCATOR MATRIX



S-LMht

S-LRI

ZYGOMA DENTAL IMPLANT SYSTEM

ZYGOMA COUPLE ABUTMENT



3 mm	S-ZCA4330
5 mm	S-ZCA4350

System » Zygoma
Material » Ti Grade 5 ELI
Torque Value » 25 Ncm
Driver » Hex Screw Driver

ANGLED ABUTMENT



3 mm	S-ZAA43203
5 mm	S-ZAA43205

System » Zygoma
Material » Ti Grade 5 ELI
Torque Value » 25 Ncm
Driver » Hex Screw Driver

MULTI UNIT ABUTMENT



ø4.3 (17°)

0.5 mm	S-ZMUA431705
1.0 mm	S-ZMUA43171
2.0 mm	S-ZMUA43172
3.0 mm	S-ZMUA43173
5.0 mm	S-ZMUA43175

System » Zygoma
Material » Ti Grade 5 ELI
Torque Value » 25 Ncm
Driver » Hex Screw Driver

LAB ANALOG



S-ZLA

System » Zygoma // Material » Ti Grade 5 ELI

HEALING CAP



3 mm	S-ZHA4330
5 mm	S-ZHA4350

System » Zygoma // Torque Value » 25 Ncm

IMPRESSION COPING PICK UP



S-ZICP

System » Zygoma // Material » Ti Grade 5 ELI

PLASTIC CYLINDER



S-ZPICYN

System » Zygoma // Material » POM

» **DIGITAL
SOLUTIONS**



BONE LEVEL

SCAN BODY



only for ø3.3 implant
S-PMSB



S-PSB

System » Digital System
Driver » Hex Screw Driver

DIGITAL ANALOG



only for ø3.3 implant
S-BMFDLA



S-BFDLA

System » Digital System
Driver » Hex Screw Driver

PREMILL



only for ø3.3 implant
ø10 mm S-BMPA10
ø14 mm S-BMPA14
ø16 mm -



S-BPA10
S-BPA14
S-BPA16

System » Digital System
Material » Ti Grade 5 ELI
Driver » Hex Screw Driver

TISSUE LEVEL

SCAN BODY



S-TSB



only for ø5.5 implant
S-TWSB

System » Digital System
Driver » Hex Screw Driver

DIGITAL ANALOG



S-TFDLA



only for ø5.5 implant
S-TWFDLA

System » Digital System
Driver » Hex Screw Driver

PREMILL



ø10 mm S-TPA10
ø14 mm S-TPA14
ø16 mm S-TPA16



only for ø5.5 implant
S-TWPA10
S-TWPA14
S-TWPA16

System » Digital System
Material » Ti Grade 5 ELI
Driver » Hex Screw Driver

MULTI UNIT ABUTMENT

SCAN BODY



S-PMUSB

System » Digital System
Driver » Hex Screw Driver

DIGITAL ANALOG



S-BMDLA

System » Digital System
Driver » Hex Screw Driver

ZR BASE



S-MUTB

System » Digital System
Driver » Hex Screw Driver

OCTA ABUTMENT

SCAN BODY



S-POSB

System » Digital System
Driver » Hex Screw Driver

DIGITAL ANALOG



S-BODLA

System » Digital System
Driver » Hex Screw Driver

ZR BASE



S-OTB

System » Digital System
Driver » Hex Screw Driver

Ti-BASE ABUTMENT

Ti-BASE SCAN BODY



only for ø3.3 implant
S-TIMSB



S-TISB

System » Digital System
Driver » Hex Screw Driver

Ti-BASE ABUTMENT



only for ø3.3 implant
S-BMTIA05



S-BTIA05

System » Digital System
Driver » Hex Screw Driver

INDEX CODE

PRODUCT	PLATFORM	HEIGHT	CODE
---------	----------	--------	------

BONE LEVEL FIXTURES

BONE LEVEL MINI FIXTURE

Bone Level Mini Fixture Straight SRA	ø 3.3	8 mm	S-BMFSR3308
Bone Level Mini Fixture Straight SRA	ø 3.3	10 mm	S-BMFSR3310
Bone Level Mini Fixture Straight SRA	ø 3.3	12 mm	S-BMFSR3312
Bone Level Mini Fixture Straight SRA	ø 3.3	14 mm	S-BMFSR3314

BONE LEVEL SHORT FIXTURE

Bone Level Short Fixture Straight SRA	ø 4.3	5 mm	S-BMFSR4305
Bone Level Short Fixture Straight SRA	ø 4.3	6 mm	S-BMFSR4306

Bone Level Short Fixture Straight SRA	ø 4.8	5 mm	S-BMFSR4805
Bone Level Short Fixture Straight SRA	ø 4.8	6 mm	S-BMFSR4806

Bone Level Short Fixture Straight SRA	ø 5.5	5 mm	S-BMFSR5505
Bone Level Short Fixture Straight SRA	ø 5.5	6 mm	S-BMFSR5506

BONE LEVEL FIXTURE

Bone Level Wide Fixture Hybrid SRA	ø 3.7	8 mm	S-BFHR3708
Bone Level Wide Fixture Hybrid SRA	ø 3.7	10 mm	S-BFHR3710
Bone Level Wide Fixture Hybrid SRA	ø 3.7	12 mm	S-BFHR3712
Bone Level Wide Fixture Hybrid SRA	ø 3.7	14 mm	S-BFHR3714

Bone Level Wide Fixture Hybrid SRA	ø 4.3	8 mm	S-BFHR4308
Bone Level Wide Fixture Hybrid SRA	ø 4.3	10 mm	S-BFHR4310
Bone Level Wide Fixture Hybrid SRA	ø 4.3	12 mm	S-BFHR4312
Bone Level Wide Fixture Hybrid SRA	ø 4.3	14 mm	S-BFHR4314

Bone Level Wide Fixture Hybrid SRA	ø 4.8	8 mm	S-BFHR4808
Bone Level Wide Fixture Hybrid SRA	ø 4.8	10 mm	S-BFHR4810
Bone Level Wide Fixture Hybrid SRA	ø 4.8	12 mm	S-BFHR4812
Bone Level Wide Fixture Hybrid SRA	ø 4.8	14 mm	S-BFHR4814

BONE LEVEL WIDE FIXTURE

Bone Level Wide Fixture Hybrid SRA	ø 5.5	8 mm	S-BWFHR5508
Bone Level Wide Fixture Hybrid SRA	ø 5.5	10 mm	S-BWFHR5510
Bone Level Wide Fixture Hybrid SRA	ø 5.5	12 mm	S-BWFHR5512

HEALING ABUTMENTS

Healing Abutment	ø 3.7	2.5 mm	S-BMHA3725
Healing Abutment	ø 3.7	3.5 mm	S-BMHA3735
Healing Abutment	ø 3.7	5 mm	S-BMHA3750
Healing Abutment	ø 3.7	7 mm	S-BMHA3770

Healing Abutment	ø 4.5	1 mm	S-BHA4510
Healing Abutment	ø 4.5	2 mm	S-BHA4520

PRODUCT	PLATFORM	HEIGHT	CODE
---------	----------	--------	------

Healing Abutment	ø 4.5	3 mm	S-BHA4530
Healing Abutment	ø 4.5	4 mm	S-BHA4540
Healing Abutment	ø 4.5	5 mm	S-BHA4550

Healing Abutment	ø 5.5	1 mm	S-BHA5510
Healing Abutment	ø 5.5	2 mm	S-BHA5520
Healing Abutment	ø 5.5	3 mm	S-BHA5530
Healing Abutment	ø 5.5	4 mm	S-BHA5540
Healing Abutment	ø 5.5	5 mm	S-BHA5550

Healing Abutment	ø 6.5	1 mm	S-BHA6510
Healing Abutment	ø 6.5	2 mm	S-BHA6520
Healing Abutment	ø 6.5	3 mm	S-BHA6530
Healing Abutment	ø 6.5	4 mm	S-BHA6540
Healing Abutment	ø 6.5	5 mm	S-BHA6550

Mini Healing Abutment for Short Implant	ø 5.0	2 mm	S-BMHA5020
Mini Healing Abutment for Short Implant	ø 5.0	3 mm	S-BMHA5030

Mini Healing Abutment for Short Implant	ø 6.0	2 mm	S-BMHA6020
Mini Healing Abutment for Short Implant	ø 6.0	3 mm	S-BMHA6030

Mini Healing Abutment for Short Implant	ø 8.0	2 mm	S-BMHA8020
Mini Healing Abutment for Short Implant	ø 8.0	3 mm	S-BMHA8030

BONE LEVEL COUPLE SYSTEM

COUPLE ABUTMENTS

Mini Couple Abutment H7	ø 3.7	1 mm	S-BMCA37107
Mini Couple Abutment H7	ø 3.7	2 mm	S-BMCA37207
Mini Couple Abutment H7	ø 3.7	3 mm	S-BMCA37307
Mini Couple Abutment H7	ø 3.7	4 mm	S-BMCA37407
Mini Couple Abutment H7	ø 3.7	5 mm	S-BMCA37507

Couple Abutment H7	ø 4.5	1 mm	S-BCA45107
Couple Abutment H7	ø 4.5	2 mm	S-BCA45207
Couple Abutment H7	ø 4.5	3 mm	S-BCA45307
Couple Abutment H7	ø 4.5	4 mm	S-BCA45407
Couple Abutment H7	ø 4.5	5 mm	S-BCA45507

Couple Abutment H7	ø 5.5	1 mm	S-BCA55107
Couple Abutment H7	ø 5.5	2 mm	S-BCA55207
Couple Abutment H7	ø 5.5	3 mm	S-BCA55307
Couple Abutment H7	ø 5.5	4 mm	S-BCA55407
Couple Abutment H7	ø 5.5	5 mm	S-BCA55507

Couple Abutment H7	ø 6.5	1 mm	S-BCA65107
Couple Abutment H7	ø 6.5	2 mm	S-BCA65207
Couple Abutment H7	ø 6.5	3 mm	S-BCA65307
Couple Abutment H7	ø 6.5	4 mm	S-BCA65407
Couple Abutment H7	ø 6.5	5 mm	S-BCA65507

PRODUCT	PLATFORM	HEIGHT	CODE
Mini Couple Abutment for Short Implant	∅ 5.0	2 mm	S-BMCA5020
Mini Couple Abutment for Short Implant	∅ 5.0	3 mm	S-BMCA5030
Mini Couple Abutment for Short Implant	∅ 6.0	2 mm	S-BMCA6020
Mini Couple Abutment for Short Implant	∅ 6.0	3 mm	S-BMCA6030
Mini Couple Abutment for Short Implant	∅ 8.0	2 mm	S-BMCA8020
Mini Couple Abutment for Short Implant	∅ 8.0	3 mm	S-BMCA8030

ANGLED ABUTMENTS

Mini Angled Abutment	∅ 3.7	1 mm 15°	S-BMAA37151
Mini Angled Abutment	∅ 3.7	2 mm 15°	S-BMAA37152
Mini Angled Abutment	∅ 3.7	4 mm 15°	S-BMAA37154

Angled Abutment	∅ 4.5	2 mm 15°	S-BAA45152
Angled Abutment	∅ 4.5	2 mm 25°	S-BAA45252
Angled Abutment	∅ 4.5	4 mm 15°	S-BAA45154
Angled Abutment	∅ 4.5	4 mm 25°	S-BAA45254

Angled Abutment	∅ 5.5	2 mm 15°	S-BAA55152
Angled Abutment	∅ 5.5	2 mm 25°	S-BAA55252
Angled Abutment	∅ 5.5	4 mm 15°	S-BAA55154
Angled Abutment	∅ 5.5	4 mm 25°	S-BAA55254

Angled Abutment	∅ 6.5	2 mm 15°	S-BAA65152
Angled Abutment	∅ 6.5	2 mm 25°	S-BAA65252
Angled Abutment	∅ 6.5	4 mm 15°	S-BAA65154
Angled Abutment	∅ 6.5	4 mm 25°	S-BAA65254

TEMPORARY ABUTMENTS

Mini Temporary Abutment	∅ 3.7		S-BMTA37
Temporary Abutment	∅ 4.5		S-BTA45
Temporary Abutment	∅ 5.5		S-BTA55

IMPRESSION COPING PICK UP

Mini Impression Coping Pick up	∅ 3.7		S-BMICP37L
Mini Impression Coping Pick up	∅ 3.7		S-BMICP37S

Impression Coping Pick up	∅ 4.5		S-BICP45L
Impression Coping Pick up	∅ 4.5		S-BICP45S

Impression Coping Pick up	∅ 5.5		S-BICP55L
Impression Coping Pick up	∅ 5.5		S-BICP55S

Impression Coping Pick up	∅ 6.5		S-BICP65L
Impression Coping Pick up	∅ 6.5		S-BICP65S

Impression Coping Pick up for short	∅ 5.0		S-BMSICP50
Impression Coping Pick up for short	∅ 6.0		S-BMSICP60

PRODUCT	PLATFORM	HEIGHT	CODE
IMPRESSION COPING TRANSFER			
Mini Impression Coping Transfer	∅ 3.7		S-BMICT37L
Mini Impression Coping Transfer	∅ 3.7		S-BMICT37S
Mini Impression Coping Transfer for short	∅ 5.0		S-BMSICT50
Mini Impression Coping Transfer for short	∅ 6.0		S-BMSICT60
Impression Coping Transfer	∅ 4.5		S-BIP45

FIXTURE LAB ANALOG

Mini Fixture Lab Analog	∅ 3.3	12 mm	S-BMFLA
Fixture Lab Analog	∅ 4.3	12 mm	S-BFLA

BONE LEVEL SOLID SYSTEM

SOLID ABUTMENTS

Mini Solid Abutment	∅ 3.7	1 mm	S-BMSA3710
Mini Solid Abutment	∅ 3.7	2 mm	S-BMSA3720
Mini Solid Abutment	∅ 3.7	3 mm	S-BMSA3730
Mini Solid Abutment	∅ 3.7	4 mm	S-BMSA3740
Mini Solid Abutment	∅ 3.7	5 mm	S-BMSA3750

Solid Abutment	∅ 4.5	1 mm	S-BSA4510
Solid Abutment	∅ 4.5	2 mm	S-BSA4520
Solid Abutment	∅ 4.5	3 mm	S-BSA4530
Solid Abutment	∅ 4.5	4 mm	S-BSA4540
Solid Abutment	∅ 4.5	5 mm	S-BSA4550

Solid Abutment	∅ 5.5	1 mm	S-BSA5510
Solid Abutment	∅ 5.5	2 mm	S-BSA5520
Solid Abutment	∅ 5.5	3 mm	S-BSA5530
Solid Abutment	∅ 5.5	4 mm	S-BSA5540
Solid Abutment	∅ 5.5	5 mm	S-BSA5550

SOLID ABUTMENT IMPRESSION CAPS

Mini Solid Abutment Impression Cap	∅ 3.7		S-BMSAICap37
Solid Abutment Impression Cap	∅ 4.5		S-BSAICap45
Solid Abutment Impression Cap	∅ 5.5		S-BSAICap55

SOLID LAB ANALOG

Mini Solid Lab Analog	∅ 3.7		S-BMSALA37
Solid Lab Analog	∅ 4.5		S-BSALA45
Solid Lab Analog	∅ 5.5		S-BSALA55

BONE LEVEL OCTA SYSTEM

OCTA ABUTMENTS

Mini Octa Abutment	∅ 4.8	1 mm	S-BMOA10
Mini Octa Abutment	∅ 4.8	2 mm	S-BMOA20
Mini Octa Abutment	∅ 4.8	3 mm	S-BMOA30
Mini Octa Abutment	∅ 4.8	4 mm	S-BMOA40
Mini Octa Abutment	∅ 4.8	5 mm	S-BMOA50

PRODUCT	PLATFORM	HEIGHT	CODE
Octa Abutment	∅ 4.8	1 mm	S-BOA10
Octa Abutment	∅ 4.8	2 mm	S-BOA20
Octa Abutment	∅ 4.8	3 mm	S-BOA30
Octa Abutment	∅ 4.8	4 mm	S-BOA40
Octa Abutment	∅ 4.8	5 mm	S-BOA50

OCTA ABUTMENT HEALING CAP

Octa Abutment Healing Cap	∅ 4.9	4 mm	S-OAHCap
Octa Abutment Healing Cap Large	∅ 4.9	6 mm	S-OAHCAP-L

OCTA SYSTEM PLASTIC CYLINDER OCTA

Octa System Plastic Cylinder Octa	∅ 5.2 mm		S-OAPICy
-----------------------------------	----------	--	----------

OCTA SYSTEM PLASTIC CYLINDER NON-OCTA

Octa System Plastic Cylinder Non-Octa	∅ 5.2 mm		S-OAPICyN
---------------------------------------	----------	--	-----------

OCTA SYSTEM TITANIUM CYLINDER

Octa System Plastic Cylinder	∅ 5.2 mm		S-OATICy
------------------------------	----------	--	----------

OCTA SYSTEM TITANIUM CYLINDER NON-OCTA

Octa System Plastic Cylinder Non-Hex	∅ 5.2		S-OATICyN
--------------------------------------	-------	--	-----------

OCTA SYSTEM IMPRESSION COPING PICK UP

Octa System Impression Coping Pick Up	∅ 5.5	10 mm	S-ICOAP
---------------------------------------	-------	-------	---------

OCTA SYSTEM IMPRESSION COPING PICK UP NON-OCTA

Octa System Impression Coping Pick Up Non-Hex	∅ 5.5	10 mm	S-ICOAPN
---	-------	-------	----------

OCTA SYSTEM IMPRESSION COPING TRANSFER

Octa System Impression Coping Transfer	∅ 5.0	7 mm	S-ICOAT
--	-------	------	---------

OCTA SYSTEM IMPRESSION COPING TRANSFER NON-OCTA

Octa System Impression Coping Transfer Non-Hex	∅ 5.0	7 mm	S-ICOATN
--	-------	------	----------

OCTA ABUTMENT LAB ANALOG

Octa Abutment Lab Analog	∅ 4.8	14 mm	S-OALA
--------------------------	-------	-------	--------

BONE LEVEL O-RING SYSTEM

O-RING ABUTMENTS

Mini O-Ring Abutment	∅ 2.9	0.5 mm	S-BMORA05
Mini O-Ring Abutment	∅ 3.7	2 mm	S-BMORA20
Mini O-Ring Abutment	∅ 3.7	4 mm	S-BMORA40

O-Ring Abutment	∅ 3.42	0.5 mm	S-BORA05
O-Ring Abutment	∅ 3.42	1 mm	S-BORA10

O-Ring Abutment	∅ 4.5	2 mm	S-BORA20
O-Ring Abutment	∅ 4.5	3 mm	S-BORA30
O-Ring Abutment	∅ 4.5	4 mm	S-BORA40
O-Ring Abutment	∅ 4.5	6 mm	S-BORA60

PRODUCT	PLATFORM	HEIGHT	CODE
DIRECTIONAL RINGS			
Directional Rings			S-ORRa

RETAINERS

Retainer	∅ 5.5	2 mm	S-ORR
Retainer	∅ 5.5	3.5 mm	S-ORR-C
Retainer Pink Cap			S-ORG-PCap
Retainer Yellow Cap			S-ORG-YCap
Retainer Green Cap			S-ORG-GCap

O-RING LAB ANALOG

O-Ring Lab Analog	∅ 4.1		S-ORLA
-------------------	-------	--	--------

BONE LEVEL LOCATOR SYSTEM

LOCATOR ABUTMENTS

Locator Mini Abutment	∅ 3.86mm	1.5mm	S-BMLA15
Locator Mini Abutment	∅ 3.86mm	2.5mm	S-BMLA25
Locator Mini Abutment	∅ 3.86mm	3.5mm	S-BMLA35
Locator Mini Abutment	∅ 3.86mm	4.5mm	S-BMLA45
Locator Mini Abutment	∅ 3.86mm	5.5mm	S-BMLA55
Locator Mini Abutment	∅ 3.86mm	6.5mm	S-BMLA65

Locator Abutment	∅ 3.86mm	1.5mm	S-BLA15
Locator Abutment	∅ 3.86mm	2.5mm	S-BLA25
Locator Abutment	∅ 3.86mm	3.5mm	S-BLA35
Locator Abutment	∅ 3.86mm	4.5mm	S-BLA45
Locator Abutment	∅ 3.86mm	5.5mm	S-BLA55
Locator Abutment	∅ 3.86mm	6.5mm	S-BLA65

LOCATOR MATRIX

Locator Matrix (includes processing cap)	∅ 5.5mm		S-LMHt
Locator Cap (includes 4 pcs caps)	∅ 4.86 mm		S-LRI

LOCATOR ANALOG

Locator Lab Analog	∅ 4 mm		S-BLLA
--------------------	--------	--	--------

BONE LEVEL MULTI-UNIT SYSTEM

MULTI-UNIT ANGLED ABUTMENT

Bone Level Multi Unit Angled Abutment	∅ 4.8x17	1.0 mm	S-BMUA48171
Bone Level Multi Unit Angled Abutment	∅ 4.8x17	2.5 mm	S-BMUA48172
Bone Level Multi Unit Angled Abutment	∅ 4.8x17	3.5 mm	S-BMUA48173
Bone Level Multi Unit Angled Abutment	∅ 4.8x35	1.0 mm	S-BMUA48351
Bone Level Multi Unit Angled Abutment	∅ 4.8x35	2.5 mm	S-BMUA48352
Bone Level Multi Unit Angled Abutment	∅ 4.8x35	3.5 mm	S-BMUA48353

MULTI-UNIT ABUTMENT

Bone Level Multi Unit Mini Abutment	∅ 4.8	0.5 mm	S-BMMUA4805
Bone Level Multi Unit Mini Abutment	∅ 4.8	1.5 mm	S-BMMUA4815
Bone Level Multi Unit Mini Abutment	∅ 4.8	2.5 mm	S-BMMUA4825

PRODUCT	PLATFORM	HEIGHT	CODE
Bone Level Multi Unit Mini Abutment	ø 4.8	3.5 mm	S-BMMUA4835
Bone Level Multi Unit Mini Abutment	ø 4.8	4.5 mm	S-BMMUA4845
Bone Level Multi Unit Abutment	ø 4.8	0.5 mm	S-BMUA4805
Bone Level Multi Unit Abutment	ø 4.8	1.5 mm	S-BMUA4815
Bone Level Multi Unit Abutment	ø 4.8	2.5 mm	S-BMUA4825
Bone Level Multi Unit Abutment	ø 4.8	3.5 mm	S-BMUA4835
Bone Level Multi Unit Abutment	ø 4.8	4.5 mm	S-BMUA4845
MULTI UNIT HEALING CAP			
Multi Unit Healing Cap	ø 4.8	4.0 mm	S-MUAHCAP
Multi Unit Healing Cap	ø 4.8	6.0 mm	S-MUAHCAPL
MULTI UNIT IMPRESSION COPING PICK UP			
Bone Level Multi Unit Impression Coping Pick Up	ø 4.8		S-BMUAICP
MULTI UNIT ABUTMENT PLASTIC CYLINDER			
Multi Unit Abutment Plastic Cylinder	ø 4.8		S-MUAPICyN
MULTI UNIT ABUTMENT TITANIUM CYLINDER			
Multi Unit Abutment Titanium Cylinder	ø 4.8		S-MUATICyN
MULTI UNIT LAB ANALOG			
Multi Unit Abutment Lab Analog	ø 4.8	14 mm	S-MUALA
BAR SYSTEM			
Bar burnout			S-IB
Bar Positioning			S-IBP
Bar Housing			S-IBH
Bar Clip Yellow			S-IBClip-y
Bar Clip Pink			S-IBClip-P
Bar Clip Red			S-IBClip-r
Tool for Inserting Clip			S-BMIA
Key for Parallelometer			S-TFIC
TISSUE LEVEL FIXTURES			
TISSUE LEVEL FIXTURE HYBRID			
Tissue Level Wide Fixture Hybrid SRA	ø 3.7	8 mm	S-TFHR3708
Tissue Level Wide Fixture Hybrid SRA	ø 3.7	10 mm	S-TFHR3710
Tissue Level Wide Fixture Hybrid SRA	ø 3.7	12 mm	S-TFHR3712
Tissue Level Wide Fixture Hybrid SRA	ø 3.7	14 mm	S-TFHR3714
Tissue Level Wide Fixture Hybrid SRA	ø 4.3	8 mm	S-TFHR4308
Tissue Level Wide Fixture Hybrid SRA	ø 4.3	10 mm	S-TFHR4310
Tissue Level Wide Fixture Hybrid SRA	ø 4.3	12 mm	S-TFHR4312
Tissue Level Wide Fixture Hybrid SRA	ø 4.3	14 mm	S-TFHR4314
Tissue Level Wide Fixture Hybrid SRA	ø 4.8	8 mm	S-TFHR4808
Tissue Level Wide Fixture Hybrid SRA	ø 4.8	10 mm	S-TFHR4810

PRODUCT	PLATFORM	HEIGHT	CODE
Tissue Level Wide Fixture Hybrid SRA	ø 4.8	12 mm	S-TFHR4812
Tissue Level Wide Fixture Hybrid SRA	ø 4.8	14 mm	S-TFHR4814
TISSUE LEVEL SHORT FIXTURE STRAIGHT			
Tissue Level Short Fixture Straight SRA	ø 4.3	4 mm	S-TFSR4304
Tissue Level Short Fixture Straight SRA	ø 4.3	5 mm	S-TFSR4305
Tissue Level Short Fixture Straight SRA	ø 4.3	6 mm	S-TFSR4306
Tissue Level Short Fixture Straight SRA	ø 4.8	4 mm	S-TFSR4804
Tissue Level Short Fixture Straight SRA	ø 4.8	5 mm	S-TFSR4805
Tissue Level Short Fixture Straight SRA	ø 4.8	6 mm	S-TFSR4806
Tissue Level Short Fixture Straight SRA	ø 5.5	4 mm	S-TWFSR5504
Tissue Level Short Fixture Straight SRA	ø 5.5	5 mm	S-TWFSR5505
Tissue Level Short Fixture Straight SRA	ø 5.5	6 mm	S-TWFSR5506
TISSUE LEVEL WIDE FIXTURE HYBRID			
Tissue Level Wide Fixture Hybrid SRA	ø 5.5	8 mm	S-TWFHR5508
Tissue Level Wide Fixture Hybrid SRA	ø 5.5	10 mm	S-TWFHR5510
Tissue Level Wide Fixture Hybrid SRA	ø 5.5	12 mm	S-TWFHR5512
TISSUE LEVEL HEALING ABUTMENTS			
Tissue Level Healing Abutment	ø 5.2	3.5 mm	S-THA35
Tissue Level Healing Abutment	ø 5.2	5.5 mm	S-THA55
Tissue Level Healing Abutment	ø 6.4	3.5 mm	S-THAW35
Tissue Level Healing Abutment	ø 6.4	5.5 mm	S-THAW55
TISSUE LEVEL COUPLE SYSTEM			
TISSUE LEVEL COUPLE ABUTMENTS			
Tissue Level Couple Abutment	ø 5.2	1 mm	S-TCA10
Tissue Level Couple Abutment	ø 5.2	2 mm	S-TCA20
Tissue Level Couple Abutment	ø 5.2	3 mm	S-TCA30
Tissue Level Couple Abutment	ø 5.2	4 mm	S-TCA40
Tissue Level Couple Abutment Wide	ø 6.4	1 mm	S-TCAW10
Tissue Level Couple Abutment Wide	ø 6.4	2 mm	S-TCAW20
Tissue Level Couple Abutment Wide	ø 6.4	3 mm	S-TCAW30
Tissue Level Couple Abutment Wide	ø 6.4	4 mm	S-TCAW40
TISSUE LEVEL ANGLED ABUTMENTS			
Tissue Level Angled Abutment	ø 3.7	15°	S-TAA15
Tissue Level Angled Abutment	ø 3.7	25°	S-TAA25
TISSUE LEVEL SYNOCTA ABUTMENT OCTA			
Tissue Level Synocta Abutment Octa	ø 3.5	6 mm	S-TSOA
TISSUE LEVEL IMPRESSION COPING PICK UP OCTA			
Tissue Level Impression Coping Pick Up Octa	ø 5.5	11.75 mm	S-TICP

PRODUCT	PLATFORM	HEIGHT	CODE
Tissue Level Impression Coping Pick Up Wide Octa	ø 6.4	11.75 mm	S-TICPW
TISSUE LEVEL IMPRESSION COPING TRANSFER			
Tissue Level Impression Coping Transfer Octa	ø 5.5	9 mm	S-TICT
TISSUE LEVEL FIXTURE LAB ANALOG			
Tissue Level Fixture Lab Analog	ø 4.8	14 mm	S-TFLA
Tissue Level Fixture Large Lab Analog	ø 6.0	14 mm	S-TFLLA
TISSUE LEVEL SOLID SYSTEM			
TISSUE LEVEL SOLID ABUTMENTS			
Tissue Level Solid Abutment	ø 3.5	4 mm	S-TSA40
Tissue Level Solid Abutment	ø 3.5	5.5 mm	S-TSA55
Tissue Level Solid Abutment	ø 3.5	7 mm	S-TSA70
Tissue Level Solid Abutment Large	ø 4.5	4 mm	S-TSAL40
Tissue Level Solid Abutment Large	ø 4.5	5 mm	S-TSAL50
Tissue Level Solid Abutment Large	ø 4.5	7 mm	S-TSAL70
TISSUE LEVEL SOLID ABUTMENT POSITIONING CYLINDER			
Tissue Level Solid Abutment Positioning Cylinder	ø 10.2		S-TSAPoCy
TISSUE LEVEL SOLID ABUTMENT PLASTIC COPING			
Tissue Level Solid Plastic Coping Crown	ø 5.3	10 mm	S-TSAPICc
Tissue Level Solid Plastic Coping Bridge	ø 5.3	10 mm	S-TSAPICb
TISSUE LEVEL SOLID ABUTMENT IMPRESSION CAP			
Tissue Level Solid Abutment Impression Cap	ø 10.2		S-TSAICap
TISSUE LEVEL SOLID ABUTMENT LAB ANALOG			
Tissue Level Solid Abutment Lab Analog	ø 4.8	14 mm	S-TSALA40
Tissue Level Solid Abutment Lab Analog	ø 4.8	15.5 mm	S-TSALA55
Tissue Level Solid Abutment Lab Analog	ø 4.8	17 mm	S-TSALA70
Tissue Level Solid Abutment Large Lab Analog	ø 6.0	14 mm	S-TSALLA40
Tissue Level Solid Abutment Large Lab Analog	ø 6.0	15.5 mm	S-TSALLA55
Tissue Level Solid Abutment Large Lab Analog	ø 6.0	17 mm	S-TSALLA70
TISSUE LEVEL OCTA SYSTEM			
TISSUE LEVEL OCTA ABUTMENT			
Tissue Level Octa Abutment	ø 3.5		S-TOA
OCTA ABUTMENT HEALING CAP			
Octa Abutment Healing Cap	ø 4.9	4 mm	S-OAHCAP
Octa Abutment Healing Cap Large	ø 4.9	6 mm	S-OAHCAP-L
TISSUE LEVEL OCTA SYSTEM PLASTIC CYLINDER			
Octa System Plastic Cylinder Octa	ø 5		S-OAPICy

PRODUCT	PLATFORM	HEIGHT	CODE
TISSUE LEVEL OCTA ABUTMENT PLASTIC CYLINDER NON-OCTA			
Tissue Level Octa Abutment Plastic Cylinder Non-Octa	ø 5		S-OAPICyN
TISSUE LEVEL OCTA SYSTEM IMPRESSION COPING PICK UP			
Tissue Level Octa System Impression Coping Pick Up Octa	ø 5.5		S-ICOAP
Tissue Level Octa System Impression Coping Pick Up Non-Octa	ø 5.5		S-ICOAPN
TISSUE LEVEL OCTA SYSTEM IMPRESSION COPING TRANSFER			
Tissue Level Octa System Impression Coping Transfer Octa	ø 5.0		S-ICOAT
Tissue Level Octa System Impression Coping Transfer Non-Octa	ø 5.0		S-ICOATN
TISSUE LEVEL OCTA ABUTMENT LAB ANALOG			
Tissue Level Octa Abutment Lab Analog	ø 4.8		S-OALA
TISSUE LEVEL O-RING SYSTEM			
TISSUE LEVEL O-RING ABUTMENTS			
Tissue Level O-Ring Abutment	ø 3.5	0.5 mm	S-TORA05
Tissue Level O-Ring Abutment	ø 3.5	2 mm	S-TORA20
Tissue Level O-Ring Abutment	ø 3.5	4 mm	S-TORA40
DIRECTIONAL RINGS			
Directional Rings			S-ORRa
RETAINERS			
Retainer	ø 5.5	2 mm	S-ORR
Retainer	ø 5.5	3.5 mm	S-ORR-C
Retainer Pink Cap			S-ORG-PCap
Retainer Yellow Cap			S-ORG-YCap
Retainer Green Cap			S-ORG-GCap
TISSUE LEVEL O-RING LAB ANALOG			
Tissue Level O-Ring Lab Analog	ø 4.1		S-ORLA
TISSUE LEVEL LOCATOR SYSTEM			
TISSUE LEVEL LOCATOR ABUTMENTS			
Tissue Level Locator Abutment	ø 4.9	1.5 mm	S-TLA15
Tissue Level Locator Abutment	ø 4.9	2.5 mm	S-TLA25
Tissue Level Locator Abutment	ø 4.9	3.5 mm	S-TLA35
Tissue Level Locator Abutment	ø 4.9	4.5 mm	S-TLA45
Tissue Level Locator Abutment	ø 4.9	5.5 mm	S-TLA55
Tissue Level Locator Abutment	ø 4.9	6.5 mm	S-TLA65
LOCATOR MATRIX			
Locator Matrix (includes processing cap)	ø 5.5mm		S-LMHt
Locator Cap (includes 4 pcs caps)	ø 4.86 mm		S-LRI
TISSUE LEVEL LOCATOR LAB ANALOG			
Tissue Level Locator Lab Analog			S-TLLA

PRODUCT	PLATFORM	HEIGHT	CODE
---------	----------	--------	------

ZYGOMA IMPLANT

ZYGOMA FIXTURE

Zygoma Fixture	∅ 4.3	35 mm	S-ZYG35
Zygoma Fixture	∅ 4.3	37.5 mm	S-ZYG37.5
Zygoma Fixture	∅ 4.3	40 mm	S-ZYG40
Zygoma Fixture	∅ 4.3	42.5 mm	S-ZYG42.5
Zygoma Fixture	∅ 4.3	45 mm	S-ZYG45
Zygoma Fixture	∅ 4.3	47.5 mm	S-ZYG47.5
Zygoma Fixture	∅ 4.3	50 mm	S-ZYG50
Zygoma Fixture	∅ 4.3	52.5 mm	S-ZYG52.5
Zygoma Fixture	∅ 4.3	55 mm	S-ZYG55
Zygoma Fixture	∅ 4.3	57.5 mm	S-ZYG57.5
Zygoma Fixture	∅ 4.3	60 mm	S-ZYG60

ATROPHIC MAXILLARY ZYGOMA FIXTURE

Atrophic Maxillary Zygoma Fixture	∅ 4.3	35 mm	S-ZAM35
Atrophic Maxillary Zygoma Fixture	∅ 4.3	37.5 mm	S-ZAM37.5
Atrophic Maxillary Zygoma Fixture	∅ 4.3	40 mm	S-ZAM40
Atrophic Maxillary Zygoma Fixture	∅ 4.3	42.5 mm	S-ZAM42.5
Atrophic Maxillary Zygoma Fixture	∅ 4.3	45 mm	S-ZAM45
Atrophic Maxillary Zygoma Fixture	∅ 4.3	47.5 mm	S-ZAM47.5
Atrophic Maxillary Zygoma Fixture	∅ 4.3	50 mm	S-ZAM50
Atrophic Maxillary Zygoma Fixture	∅ 4.3	52.5 mm	S-ZAM52.5
Atrophic Maxillary Zygoma Fixture	∅ 4.3	55 mm	S-ZAM55
Atrophic Maxillary Zygoma Fixture	∅ 4.3	57.5 mm	S-ZAM57.5
Atrophic Maxillary Zygoma Fixture	∅ 4.3	60 mm	S-ZAM60

ONCOLOGICAL ZYGOMA FIXTURE

Oncological Zygoma Fixture	∅ 4.3	25 mm	SZ-ONC25
Oncological Maxillary Zygoma Fixture	∅ 4.3	27.5 mm	SZ-ONC27.5
Oncological Maxillary Zygoma Fixture	∅ 4.3	30 mm	SZ-ONC30
Oncological Maxillary Zygoma Fixture	∅ 4.3	32.5 mm	SZ-ONC32.5

ABUTMENTS

Zygoma Abutment		3 mm	S-ZCA4330
Zygoma Abutment		5 mm	S-ZCA4350

ZYGOMA ANGLED ABUTMENTS

Zygoma Angled Abutment		3 mm	S-ZAA43203
Zygoma Angled Abutment		5 mm	S-ZAA43205

ZYGOMA MULTI UNIT ABUTMENT

Zygoma Multi Unit Abutment		0.5 mm	S-ZMUA431705
Zygoma Multi Unit Abutment		1 mm	S-ZMUA43171
Zygoma Multi Unit Abutment		2 mm	S-ZMUA43172
Zygoma Multi Unit Abutment		3 mm	S-ZMUA43173
Zygoma Multi Unit Abutment		5 mm	S-ZMUA43175

PRODUCT	PLATFORM	HEIGHT	CODE
---------	----------	--------	------

CAPS

Zygoma Healing Cap		3 mm	S-ZHA4330
Zygoma Healing Cap		5 mm	S-ZHA4350

COPING

Zygoma Impression Coping Pick Up			S-ZICP
----------------------------------	--	--	--------

PLASTIC CYLINDER

Zygoma Plastic Cylinder			S-ZPICYN
-------------------------	--	--	----------

LAB ANALOG

Zygoma Lab Analog			S-ZLA
-------------------	--	--	-------

DRIVERS

Implant Drivers (For ratchet and contra angle)

Bone Level Mini Fixture Driver Machine Short			S-BMFDMS
Bone Level Mini Fixture Driver Machine Long			S-BMFDML

Bone Level Mini Fixture Driver Ratchet Short			S-BMFDMS
Bone Level Mini Fixture Driver Ratchet Long			S-BMFDRL

Bone Level Fixture Driver Machine Short			S-BFDMS
Bone Level Fixture Driver Machine Long			S-BFDML

Bone Level Fixture Driver Ratchet Short			S-BFDMS
Bone Level Fixture Driver Ratchet Long			S-BFDRL

Tissue Level Fixture Driver Ratchet Short			S-TFDMS
Tissue Level Fixture Driver Ratchet Long			S-TFDRL

Tissue Level Fixture Driver Machine Short			S-TFDMS
Tissue Level Fixture Driver Machine Long			S-TFDML

Hex driver (extra short, short, long, extra long)			S-HDRS / S-HDRL
---	--	--	-----------------

Hex Driver Contra Angle (extra short, short, long, extra long)			S-HDRS / S-HDRL
--	--	--	-----------------

O-ring Driver			S-ORADS/L
---------------	--	--	-----------

Octa Driver			S-OADS/L
-------------	--	--	----------

Solid Driver			S-SAD
--------------	--	--	-------

ORRO Insert Adaptor			S-OIA
---------------------	--	--	-------

Bar Matrix Insert Adaptor			S-BMIA
---------------------------	--	--	--------

PRODUCT	PLATFORM	HEIGHT	CODE
SCREWS			
CLOSING SCREW			
BONE LEVEL			
Bone Level Mini Closing Screw	∅ 2.8	4.5 mm	S-BMCS
Bone Level Closing Screw	∅ 3.36	6.5 mm	S-BCS
TISSUE LEVEL			
Tissue Level Closing-Screw Small	∅ 3.5	5.35 mm	S-TCSS
Tissue Level Closing-Screw Large	∅ 5.2	6.4mm	S-TCSL
Tissue Level Closing-Screw Wide	∅ 6.4	6.4 mm	S-TCSW
ABUTMENT SCREWS			
BONE LEVEL			
Bone Level Mini Abutment Screw	∅ 2.1	7.64 mm	S-BMAS
Bone Level Abutment Screw	∅ 2.3	9.32 mm	S-BAS
Bone Level Multi-Unit Abutment Screw	∅ 1.96	9 mm	S-BMUAS
TISSUE LEVEL			
Tissue Level Abutment Screw	∅ 2.5	7.1	S-TAS
COMMON SCREWS			
Octa Abutment Screw	∅ 2.5	4.5 mm	S-OAS
Multi-Unit Abutment Healing Screw	∅ 2.2	4 mm	S-MUAHS
DIGITAL SOLUTIONS			
BONE LEVEL SCAN BODY			
Bone Level Mini Scan Body			S-PMSB
Bone Level Scan Body			S-PSB
BONE LEVEL DIGITAL ANALOG			
Bone Level Mini Digital Analog			S-BMFDLA
Bone Level Digital Analog			BFDLA
BONE LEVEL PREMILL			
Bone Level Mini Premill	10 mm		S-BMPA10
Bone Level Mini Premill	14 mm		S-BMPA14
Bone Level Premill	10 mm		S-BPA10
Bone Level Premill	14 mm		S-BPA14
Bone Level Premill	16 mm		S-BPA16
TISSUE LEVEL SCAN BODY			
Tissue Level Mini Scan Body			S-TSB
Tissue Level Scan Body			S-TWFDLA

PRODUCT	PLATFORM	HEIGHT	CODE
TISSUE LEVEL DIGITAL ANALOG			
Tissue Level Mini Digital Analog			S-TFDLA
Tissue Level Digital Analog			S-TWFDLA
TISSUE LEVEL PREMILL			
Tissue Level Mini Premill	10 mm		S-TPA10
Tissue Level Mini Premill	14 mm		S-TPA14
Tissue Level Mini Premill	16 mm		S-TPA16
Tissue Level Premill	10 mm		S-TWPA10
Tissue Level Premill	14 mm		S-TWPA14
Tissue Level Premill	16 mm		S-TWPA16
MULTI UNIT ABUTMENT			
Multi Unit Abutment Scan Body			S-PMUSB
Multi Unit Abutment Digital Analog			S-BMDLA
Multi Unit Abutment ZR Base			S-MUTB
OCTA ABUTMENT			
Octa Abutment Scan Body			S-POSB
Octa Abutment Digital Analog			S-BODLA
Octa Abutment ZR Base			S-OTB
Ti-BASE ABUTMENT			
Ti-Base Mini Abutment			S-BMTIA05
Ti-Base Abutment			S-BTIA05
Ti-BASE SCAN BODY			
Ti-Base Mini Scan Body			S-TIMSB
Ti-Base Scan Body			S-TISB

